




CLEAR BLUE SMILES™

SureSmile Aligner™

Ordering Replacement Aligners

1. Login into SureSmile



The screenshot shows the SureSmile login interface. At the top left is the Dentsply Sirona logo. The header area contains the SureSmile® logo. Below the header, the text "Log in with your account details" is centered. A light gray box contains the login form with two input fields: "Email address" and "Password", followed by a "Login" button. Below the form, there are two links: "Forgot your password?" and "Looking to add a new user to your account?". To the right of these links is a blue button that says "New to SureSmile*? Click here to create an account". At the bottom left, there is a copyright notice: "© 2022 OraMetrix. All rights reserved.". At the bottom right, there is a link for "Privacy Policy".

Dentsply Sirona

SureSmile®

Log in with your account details

Email address

Password

Login

[Forgot your password?](#)


[Looking to add a new user to your account?](#)

New to SureSmile*?
Click here to create an account

© 2022 OraMetrix. All rights reserved.

[Privacy Policy](#)

2. Select on Patients tab



The screenshot displays the SureSmile Aligner software interface. At the top, there is a header bar with a logo, a 'Patient Database Tasks' section, a search bar for 'Patient ID', a 'New Patient' button, and a user profile for 'Dr. Brian Gaudreault SOUTHEAST ORTHO-ALIGNER'. Below the header, there are filters for 'Sort by' (Task), 'Show Tasks of' (All Types), and 'Tasks for' (Dr. Brian Gau...). A red arrow points to the 'Patients' tab in the navigation bar, which is currently selected. The main content area shows a table with patient and task information.

Flag	Card	Last name	First name	Patient ID	Owner	Item	Task	Due Date	Notes	Setup Approval Date	Status	Deferred To
■	👁	Hudzik	Anna	G8S7	Dr. Brian Gaudreault		Review setup & staged models ⓘ	2022-08-01	📄	2022-08-19	Active	

3. Select on Setup & Staging

The screenshot displays the 'Patients (57)' tab in the CLEAR BLUE SMILES software. The interface includes a top navigation bar with tabs for Tasks, Patients (57), Shipments, Jobs, and Reports. A 'Sort by' dropdown is set to 'Modification ...'. The 'Patients for' dropdown is set to 'All', and the 'Show' dropdown is set to 'Active Patien...'. The 'View' button is set to a grid view.

The patient grid contains 12 cards, each with a profile picture, name, ID, and status. A red arrow points to the card for 'Albers, Destinee'.

Profile Picture	Name	ID	Status
	Patient, New	ID: R3F3	Setup & Staging 1 2022-09-29 Approved / Approved
	Albers, Carrie	ID: E8J5	
	Albers, Destinee	ID: W4X4	Setup & Staging 1 2022-09-29 Approved / Approved
	test222, test	ID: N4G6	
	Lee, Brian	ID: S6X2	Setup & Staging 1 2021-08-30 Approved / Approved
	testggg, testtt	ID: T2X7	
	test, test	ID: K4H6	
	Trim, Test	ID: S7H8	
	Harris, Kevin	ID: R3E7	Setup & Staging 1 2022-09-27 Approved / Approved
	Ferguson, Carrie	ID: L6K9	
	Davis, Oliver	ID: K8J7	
	Patient, New	ID: Y9F9	

4. Select Order Aligners

The screenshot displays the SureSmile Aligner software interface. At the top, the header includes the 'Patient Database' tab, the patient name 'Anna Hudzik (G8S7)', the 'Setup & Staging 1 (Approved / Approved)' stage, and the 'SureSmile Aligner' title. A left sidebar contains a 'Checklist and Overview' panel with sections for 'Overview' (instructions to review notes and IPR), 'Setup and Staging Review' (listing Labial, Right Buccal, Left Buccal, Upper Occlusal, and Lower Occlusal views), and a bottom section with buttons for 'Order Retainers', 'Order Aligners' (highlighted with a red arrow), and 'Order Refinement'. The main workspace shows a 3D model of a patient's teeth with a teal circle highlighting the upper front teeth. Small colored squares (blue, yellow, red) are placed on specific teeth. The bottom of the interface features a 'Camera (Patient View) Navigation' bar with playback controls, a timeline from 1 to 29, and buttons for 'Legend' and 'Tooth movement'. A 'Get report' button is located in the top right corner.

5. Select Aligners

Order Aligners ✕

☐ Upper

T1123456789101112131415

1617181920212223242526272829

☐ Lower

T1123456789101112131415

1617181920212223242526272829

☐ Legend

SureSmile® Retainers can be ordered any time after placing your initial aligner order. A reminder to "Order Retainers" will appear in your task list 6 weeks prior to end of treatment if retainers have not previously been ordered.

Trim Type & Height

CloseNext

5. Select Next

Order Aligners ×

☐ Upper

T1123456789101112131415

1617181920212223242526272829

☐ Lower

T1123456789101112131415

1617181920212223242526272829

☐ Legend

Trim Type & Height

CloseNext

←

SureSmile® Retainers can be ordered any time after placing your initial aligner order. A reminder to "Order Retainers" will appear in your task list 6 weeks prior to end of treatment if retainers have not previously been ordered.

5. Select OK

