

CLEAR BLUE SMILES^M Why SureSmile[®] Aligners

The SureSmile Difference

- Accurate Digital Patient Registration
- Software Algorithm
- Aligner Material and Trim Line


¹ According to SureSmile doctors. Additional information on platform efficiency and effectiveness can be read online here: <u>Sachdeva, R.C.L. (2014). A review of suresmile: Efficiency and</u> <u>effectiveness – sprinkled with a dos of self reflection. Orthodontic Practice.</u>

CLEAR BLUE SMILES™

CLEAR BLUE SMILES™

Creating the Digital Patient


Creating the Digital Patient



Digital Model

© 2022 Clear Blue Smiles, Inc.


Not all Patients are Created Equal



Digital Model ≠ Digital Patient

© 2022 Clear Blue Smiles, Inc.


Anatomically Accurate Models




Registered by SureSmile Digital Lab

© 2022 Clear Blue Smiles, Inc.

How Does It Work?






Occlusal Plane

Midline and Smile Arc


CLEAR BLUE SMILES™

Leads to Accurate Smile Design



CLEAR BLUE SMILES™

Leads to Accurate Occlusal Plane



© 2022 Clear Blue Smiles, Inc.

Which Model is Correctly Oriented?





IOS deliver their digital model oriented to the Horizon



CLEAR BLUE SMILES™

IOS deliver their digital model oriented to the Horizon







SureSmile delivers their digital model oriented correctly!






Does it Change Your TX Plan and Smile Design?








CLEAR BLUE SMILES™

Does it Change Your TX Plan and Smile Design?



© 2022 Clear Blue Smiles, Inc.

CLEAR BLUE SMILESTM

SureSmile[®] Software

SureSmile® Aligners Software

- Software enables users to reduce the potential number of collisions and antagonistic movements
- You can set the software to stage IPR based on visit intervals
- Will lead to a predictable clinical protocol

¹ According to SureSmile doctors. Additional information on platform efficiency and effectiveness can be read online here: <u>Sachdeva, R.C.L. (2014). A review of suresmile: Efficiency and</u> <u>effectiveness – sprinkled with a dos of self reflection. Orthodontic Practice.</u>

CLEAR BLUE SMILES™

CLEAR BLUE SMILES™

SureSmile[®] Aligners and Attachments

• Stronger Plastic (1x-4x)¹

¹ According to SureSmile doctors. Additional information on platform efficiency and effectiveness can be read online here: <u>Sachdeva, R.C.L. (2014). A review of suresmile: Efficiency and</u> <u>effectiveness – sprinkled with a dos of self reflection. Orthodontic Practice.</u>

CLEAR BLUE SMILES™



Essix® PLUS™ Plastic

Primary use: Retainers & Bruxism Appliances

Redefining Retention

Engineered for performance, Essix® PLUS[™] Plastic combines the clarity your patients expect with the durability you can stand behind. Essix® PLUS[™] Plastic is effortlessly removed from models and trims with ease. Technicians and staff members will love the improved handling! When it comes to meeting the needs of practice and patients alike, Essix® PLUS[™] Plastic is the clear choice in retention.

rollover to enlarge image

click below for more product images



B43^{than leading competitors'}

2X MORE DURABLE

TX MORE STAIN-RESISTANT

1% MORE WEAR-RESISTANT than leading competitors'

Data on file.


CLEAR BLUE SMILES[™]


- Stronger Plastic (1x-4x)¹
- Straight Trimline


¹ According to SureSmile doctors. Additional information on platform efficiency and effectiveness can be read online here: <u>Sachdeva, R.C.L. (2014). A review of suresmile: Efficiency and</u> <u>effectiveness – sprinkled with a dos of self reflection. Orthodontic Practice.</u>

CLEAR BLUE SMILES[™]









- A) a scalloped margin,
- B) a straight cut at the gingival zenith, and
- C) a straight cut at 2 mm above the gingival zenith.

https://www.zenduradental.com/blogs/zenduradental/11688 821-trimming-differently-to-make-retention-2x-4x-stronger


CLEAR BLUE SMILES™


- Stronger Plastic (1x-4x)
- Straight Trimline (1mm)
- Plus, Software Algorithm=Fewer Attachments (66% less)¹

¹ According to SureSmile doctors. Additional information on platform efficiency and effectiveness can be read online here: <u>Sachdeva, R.C.L. (2014). A review of suresmile: Efficiency and effectiveness – sprinkled with a dos of self reflection. Orthodontic Practice.</u>

CLEAR BLUE SMILES™

Setup Alignment | Attachments





U:20 L:20



CLEAR BLUE SMILES[™]

- Stronger Plastic (1x-4x)
- Straight Trimline (1 mm)
- Plus, Software Algorithm=Fewer Attachments (66% less)
- Which leads to Shorter Tx times = Cost saving (25-40%)¹

¹ According to SureSmile doctors. Additional information on platform efficiency and effectiveness can be read online here: <u>Sachdeva, R.C.L. (2014). A review of suresmile: Efficiency and effectiveness – sprinkled with a dos of self reflection. Orthodontic Practice.</u>

CLEAR BLUE SMILES™


CLEAR BLUE SMILES™

CBCT

Premium Model Option for Bone and Roots via CBCT

© 2022 Clear Blue Smiles, Inc.


Advanced Diagnostic Metrics True 3D Premium Model



Root, Bone, and Gingiva Analysis


CLEAR BLUE SMILES™

Use CBCT to Dx and Tx with Roots in Bone





Roots in Bone Equals Stable and Healthy Periodontium



© 2022 Clear Blue Smiles, Inc.

A Digital Platform with an Anatomically Accurate 3D Model.

Better Smiles and Roots in Bone!

© 2022 Clear Blue Smiles, Inc.